

# **Miniature Stick-packaging**

# Pre-storage and Release of Reagents in Lab-on-a-Chip Systems

On-chip liquid reagent prestorage is often required in mobile analytical and diagnostic systems. Hahn-Schickard offers a novel approach for reagent pre-storage in tubularshaped composite foil pouches, called stick-packs (Fig. 1). The technology allows packaging a combination of liquid and dry reagents.

#### Release mechanism


An integrated frangible seal enables the pressure-controlled release of reagents and simplifies the layout of lab-on-a-chip (LOAC) systems, thereby making the package a functional microfluidic release unit. The frangible seal is adjusted to defined burst pressures ranging from 20 to 140 kPa. When inserted into a centrifugal microfluidic LOAC platform, the frangible seal delaminates by centrifugal forces at a defined burst frequency (Fig. 2, Fig. 3).

### One-/two-chamber stick-packs

Stick-packs featuring one chamber are for the pre-storage and release of liquid reagents ranging from 50-1200 µl. Since smaller amounts can't be released independently, liquids or solids ranging from 1-50 µl / mg are packaged in the second chamber of two-chamber stick-packs. This allows the release of



Fig. 1: Miniaturized stick-packs for combined storage and release of dry and liquid reagents



**Fig. 2:** Miniaturized stick-packs. (a) frangible seal before, and (b) during release of liquid reagents upon centrifugation

both reagents by delamination of the center and outer frangible seal subsequently.

## Attributes

Phone:

E-mail:

Website:

Fax:


Stick-packs have been successfully tested applying recovery tests, in which 99  $\pm 1\%$  of 250  $\mu$ l pre-stored liquid is released. Long-term storage tests indicat-

+49 761 203-73242

+49 761 203-73299

Info@Hahn-Schickard.de

www.Hahn-Schickard.de



**Fig. 3:** Burst frequencies in a centrifuge at variable radial positions. The marked areas correspond to standard deviations.

ed a loss of only < 0.5% for simulated 2 years storage. Air transport simulation tests were also conducted successfully. In conclusion, the developed technology enables the long-term storage and on-demand release of liquids or combinations of liquid and dry reagents.

Phone: +49 761 203-73205